



Combi-2

MEASUREMENT | TESTING | STIMULATION



Versatility of measurement, testing and stimulation in one device – the premium middle class for demanding users.

Flexibility and reliable performance are characteristics of the Combi-2 and indispensable for modern EAV practices. The versatile characteristics of the device offer the user a basis for a nature bioenergetic treatment plan. The Combi-2 combines modern metrological prerequisites for reliable treatment with various applications through low frequency electrical impulses for regulating the body's energetic meridians.



Combi-2

Measurement | Testing | Stimulation



- *Measurement system which meets the official standards of the International Medical Association for the Electroacupuncture acc. to Dr. Voll (EAV)
- *Large, removable display for individual placement
- Digital display (LCD) for indicator drop, conductance, frequency and tooth potential
- *Ergonomically designed stylus with mounted, pressure compensated, replaceable sensor tip
- *Optical point locating „OPTO-Point“ with 10-digit-LED-bar in the stylus and individual skin moisture compensation
- Large, integrated honeycomb
- Acoustic measurement support
- *Spread of scale (TP-Regler)
- Modern membrane keyboard
- Built-in high performance rechargeable battery for 25 hours full output operation
*Optional: power supply from mains with medical power adaptor
- Automatic-recharger or medical power adaptor
- Adaptable with software program **Homopath®** (www.homopath.de) and **Remifox**

***Kindling Special!** You will only find these features in this price range at Kindling!

Applications

- Measurements and evaluations according to the standards of EAV
- Testing with various products, incompatibilities, location contaminant loads and substances, etc.
- Electrical stimulation treatment: Fixed frequency (0.1 – 19.9 Hz), frequency variation (1 – 10.0 Hz), current types: sedating, tonifying, sweep waves
- Measure of electrical potential in teeth and comfortable tooth sensitivity test

Special Supplemental Options

- Expandable with all Kindling bioresonance devices
- Can be combined with EAV software **Homopath®** (www.homopath.de) and **Remifox**
- Quadrant switch for testing with hand, foot and forehead electrodes
- Extensive special accessories
- Optional: power supply from mains with medical power supply



Scope of Delivery

Pair of hand electrodes, brass
2 wires for hand electrodes
Stylus with 4-slice-tip
Ball tip
Honeycomb ø ca. 10 cm
Pair of dental electrodes
Charger, optional: power supply
Operating instructions

Dimensions (WxHxD):
21x23x18cm
Weight ca. 1700g

